



Process Chamber with Point Nozzle Scanning

Radian 7

Automated Photomask Surface Analysis



Radian 7 ICPMS with Photomask Load Ports

Author: Kyle Uhlmeyer¹, Sergey Mezheny²


Evaluation of Ultra-Trace Metal Detection Limits on Photomasks Using Radian 7 ICPMS

Synopsis


Photomasks play a critical role in defining circuit features in semiconductor lithography. Even single-digit ppt contamination of critical elements such as silver can cause defects. Radian 7 ICPMS enables ultra-trace metal detection on photomask blanks (raw quartz) as well as implanted and patterned photomasks through a fully automated scanning process and metals detection by ICPMS. For this work photomask blanks are automatically analyzed using a Radian 7 equipped with three SMIF load ports

with robotic reticle handling. A variety of scanning patterns and edge exclusion scans adapted to mask type provides detection of spatial distribution of metal and metal nanoparticle contamination on each photomask. Detection limits and recovery efficiency of metal contamination were determined by scanning a blank photomask multiple times. High droplet retention was achieved across all scans, allowing detection of trace metals at single-digit ppt levels, corresponding to E5-E8 atoms/cm² on the surface.


Sector Scanning with Radian 7 on Photomasks




Single Point (1, 1)



Line



Center Section



Whole Mask


Radian 7 includes a sector scanning function that allows users to define custom scan geometries such as single points, lines, boxes, or whole surface. Parameters such as edge exclusion, droplet overlap, and scan speed are software-controlled, enabling precise targeting of areas most sensitive to yield-impacting defects. This flexibility supports both routine QA and process development for advanced photomasks.



1. Elemental Scientific Inc.
2. Intel

Accurate, Efficient Recovery and Low Detection Limit Analysis

Accurate and Efficient Metals Analysis



Photomasks were spiked with known concentrations of trace metals and scanned under standard Radian VPD operating conditions. The recovered droplets were analyzed by ICPMS, with measured concentrations yielding an average recovery accuracy of 99.7% and a relative standard deviation (RSD) of 9.4%. First-scan recovery efficiency averaged 97.6% with an RSD of 3.4%, demonstrating the system’s reliability for quantitative analysis of surface contaminants.

Low Level Detection Limits

Element	On Mask Detection Limit (ppt)	On Mask Detection Limit (atoms/cm ²)	Element	On Mask Detection Limit (ppt)	On Mask Detection Limit (atoms/cm ²)
Li	0.2	2.2 E+07	Co	0.7	9.4 E+06
Mg	1.5	2.6 E+08	Sr	0.4	3.0 E+06
K	2.0	1.2 E+08	Mo	0.5	3.8 E+06
Ca	2.1	1.1 E+08	Ag	0.1	1.2 E+06
Ti	1.1	7.2 E+07	Cd	0.5	3.6 E+06
V	0.6	9.1 E+06	Ta	0.2	7.4 E+05
Cr	0.8	8.4 E+07	W	0.7	2.7 E+06
Mn	0.7	9.6 E+06	Pb	0.7	2.4 E+06

A blank quartz photomask was scanned seven times using a 500 µL droplet, where the recovered droplets were analyzed by ICPMS. Based on 3 σ statistical analysis of replicate results, the system achieved single-digit ppt detection limits, corresponding to E5-E8 atoms/cm² on the photomask surface. While these results are higher compared to bare silicon wafers due to impurities in the raw quartz, the detection limits are still sufficient to support trace element detection on photomasks.

