

## TRAXStation Filtering

Produces >1400  
autofiltrations  
(0.45  $\mu\text{m}$ ) per 24 hour  
day, alleviating lab  
staff fatigue



TRAXStation

Author: Tyler Herek


## Automating Filtration: Uncapping, Sampling and Consistent Syringe-Driven Autofiltration Using 13 mm Filters

### Synopsis

TRAXStation with LuerProbe automatically filters samples using any size luer filter. This study evaluates back pressure profiles and consistency of filtrate volume produced while auto-filtrating samples using 13 mm 0.45  $\mu\text{m}$  Luer filters. The automated high-pressure syringe-driven filtration process generates more consistent pressure than manual filtration, permitting use of 13 mm filters in place of conventional 33 mm 0.45  $\mu\text{m}$  filters. TRAXStation monitors back pressure in real time and adjusts syringe speed to compensate for higher pressures generated when samples with high levels of suspended solids are filtered. Advantages of 13 mm filters include lower cost, higher pressure rating, and capacity of up to 480 filters at a time.

### Maximum 0.45 $\mu\text{m}$ Filter Capacity

		Filter Diameter	
		33 mm	13 mm
TRAXStation Model	2x2x1	42	120
	4x2x2	84	240
	4x4x2	168	480




TRAXStation automatically and completely prepares samples for analysis. LuerProbe is used for sampling and filtering steps. Other probe options are available.




# Automatic Filtration by TRAXStation with LuerProbe

## Comparison of Pressure Stability in Manual Filtration vs. TRAXStation Syringe-driven Autofiltration Using 13 mm 0.45 µm Filters



Comparison of 100 pressure measurements (kPa) taken at regular intervals during the filtration of a single sample using 13 mm 0.45 µm filters – manual filtration versus TRAXStation. TRAXStation’s syringe-driven autofiltration generates higher and more consistent pressure, exceeding the pressure achieved through manual filtration even under maximum exertion. This results in faster sample-to-sample filtration times and demonstrates the enhanced pressure control and reliability of TRAXStation, which consistently produces the desired filtrate volume, improves throughput, and alleviates fatigue on lab staff.

## Reproducibility of Final Filtrate Volume Following Autofiltration Using 13 mm 0.45 µm Filters



Automatic filtration of 180 samples using 13 mm 0.45 µm filters automated using TRAXStation with LuerProbe. 12 mL of sample was initially taken up, then the filter was attached, and a portion dispensed to waste. The target filtrate volume in the tube was 10 mL, with an average actual volume of 10.05 mL and an RSD of 0.64%, indicating high reproducibility of the automated filtration process.



© Elemental Scientific |

sales@icpms.com |

www.icpms.com