

## TRAXStation Preserving

Produces racks of preserved and mixed samples as outlined in EPA Method 200.2 Section 8.1.




Authors: Davis Reed, Tyler Herek

## Evaluation of TRAXStation for Preserving Environmental Samples for Metals Analysis

### Synopsis

This study uses TRAXStation 442 to automatically preserve samples as outlined in EPA Method 200.2, applicable for groundwater, surface water, drinking water, wastewater and other environmental samples. The automation process includes uncapping, barcode reading, acid spiking, mixing, tightly recapping, and

reracking. Automatic syringe-driven addition of 1+1 HNO<sub>3</sub> spiked with cobalt was verified by measuring sample pH sixteen hours after preservation. Reproducibility was confirmed by monitoring the element Co by ICPMS. TRAXStation preserves more than one sample per minute or over 1500 samples per 24 hour day.




TRAXStation automatically and completely prepares racked samples for preservation. LuerProbe is used for spiking and mixing steps. Other probe options are available.




# Automatic Sample Preservation by TRAXStation

## Reproducibility of 0.15 mL Spikes of 1+1 HNO<sub>3</sub> into 120 Water Samples of 50 mL each



Automated preservation of 120 environmental samples using the TRAXStation. Each 50 mL water sample was spiked with 0.15 mL of 1+1 HNO<sub>3</sub> (~35% HNO<sub>3</sub>) containing 1000 ppb of cobalt. Cobalt was added solely to generate a signal on an ICPMS proportional to the amount of acid spike. As analyzed, the cobalt concentration was expected to be 3 ppb with recovery expressed as the measured concentration divided by the expected concentration. The average spike recovery was 100.06% with an RSD of 1.56%, confirming the reproducibility of the automated sample preservation system.

## pH of 120 Preserved Samples Measured 16 Hours after HNO<sub>3</sub> Spiking



After automated preservation using TRAXStation, 120 samples were held for 16 hours before pH measurement. These data show that the pH values for all samples were less than 2 (average pH = 1.86 with an RSD of 1.39%), proving that TRAXStation effectively preserves water samples as described in EPA Method 200.2, Section 8.1.

