

microFAST SingleCell

Complete Solution for Single Particle and Single Cell ICPMS Applications

Automated sample introduction system for NexION ICPMS




 **Learn More**

Biological Size Scale

Biological Size Scale

The microFAST SingleCell automation system delivers single nanoparticles or single cells intact to the ICPMS. Once the NPs or cells enter the plasma they are vaporized, atomized, and ionized. The resulting elemental content is detected by the mass spectrometer.



Single Cell ICPMS

The ability to introduce single cells into an ICPMS and measure the elemental content in each cell, or tagged to each cell, accurately takes a dedicated, well-designed sample introduction system. Having this ability allows for investigators to better understand how much of a specific nanoparticle, metallodrug, or metal-based compound enters the cell. These cells or nanoparticles will vary in size from a few nm's up to a few 100 μ m's. The typical cell types of interest will vary and with that the stability of the cell-line also varies, such that a gentle, controlled nebulization must be employed in order to not disrupt or lyse the cell.

Single Cell ICPMS Requirements

- Flexible sample volumes – μ L to mL of sample
- Ensure cells stay intact, no cell lysing
- Low pressure sample introduction
- High transport efficiency

Elemental Scientific has developed a complete sample introduction system designed specifically for single cell and nanoparticle applications. This system consists of:

- **microFAST SingleCell Autosampler**
- **CytoNeb** – single cell nebulizer
- **CytoSpray** – single cell spray chamber



microFAST SingleCell





microFAST SingleCell Automated
Sample Introduction System for NexION ICPMS
Part Number: MF-SC2-64

CytoNeb and CytoSpray

CytoNeb


- Meinhard high efficiency nebulizer
- Designed to efficiently nebulize single cells without cell rupturing
- Low internal volume
- Low backpressure (1-50 $\mu\text{L}/\text{min}$ = <50 psi)
- Low dead volume
- Patented inert PFA quick connects for nebulizer gas and samples lines




CytoNeb and CytoSpray

CytoSpray


- Spray chamber designed specifically for single cell and nanoparticle applications
- High-transport efficiency
- Separate make-up gas for better transport efficiency
- Includes one-piece ICPMS torch for simple and direct connection to the CytoSpray



microFAST SingleCell System




microFAST SingleCell System



The microFAST SingleCell system has been built for performance by optimizing the inner diameter (ID) and line lengths to ensure a quick sample transfer from vial to ICP torch.

- Fast sample-to-sample times. For example, at 20 $\mu\text{L}/\text{min}$ flow rate:
 - <3 min, when utilizing a 30 s ICPMS measurement time
 - <4 min, when utilizing a 100 s ICPMS measurement time
- High-flow sample loop washout
- Simple conversion for total metal analysis using FAST system
 - Vacuum or syringe sample loading
 - Micro or large sample volume capabilities



NexION Method Setup

Analysis Options

Setup: Advanced Options


NexION Method Setup

Sample List


	Sample	Analyte	Most Freq. Size (nm)	Mean Size (nm)	Median Size (nm)	% <100nm	Size LOD (nm)	No. of Peaks	Mean Inten. (counts)	Part. Conc. (parts/mL)	Diss. Inten. (counts)	Diss. Con (ppb)
<input type="checkbox"/>	STD1	Pt 195									3.34	
<input checked="" type="checkbox"/>	STD1	Pt 195						890	53.55		0	
<input type="checkbox"/>	Sample	Pt 195	30	32	31	100	15	2984	13.28		0.01	0.014
<input type="checkbox"/>	Sample	Pt 195	31	33	31	100	15	3181	13.59		0	0.012
<input type="checkbox"/>	30nm Sample-3	Pt 195	31	32	31	100	15	3232	13.18		0.01	0.013
<input type="checkbox"/>	50nm Sample-1	Pt 195	53	53	53	100	15	1078	55		0	0.004
<input type="checkbox"/>	50nm Sample-2	Pt 195	54	54	54	100	15	1120	57.59		0	0.004
<input type="checkbox"/>	50nm Sample-3	Pt 195	54	54	54	100	15	1138	56.88		0	0.002
<input type="checkbox"/>	70nm Sample-1	Pt 195	80	81	81	98	15	397	186.09		0	0
<input type="checkbox"/>	70nm Sample-2	Pt 195	80	80	80	98	15	407	181.52		0	0
<input type="checkbox"/>	70nm Sample-3	Pt 195	81	81	81	97	15	372	188.16		0	-0.001
<input type="checkbox"/>	BLK-1	Pt 195						1			0	0.002
<input type="checkbox"/>	BLK-2	Pt 195						1			0	0.001

Example Histogram


50nm Pt	Mean Size (nm)
Run 1	53
Run 2	54
Run 3	54
Average	53.7
STD Dev.	0.6
%RSD	1.1



Single Cell Introduction Kit Performance



Backpressure was recorded for each sample flow rate using DI water, 1X PBS, or 50 mM Tris buffer as the carrier solution. Larger ID tubing can be substituted to achieve lower backpressures.




Backpressure was recorded for varying nebulizer gas flow rates using 10 and 50 µL/min sample flow (DI water as the carrier solution).

Nanoparticle Performance

50 nm Au NPs


Typical Transport Efficiency for 50 nm Au NPs = ~80% or greater



Ten 50 nm Au NPs were prepared under the same conditions and analyzed over an 8-day period. Samples were sonicated before each day's analysis. Data points represent the average response for each sample over the 8 days. Error bars represent ± 1 standard deviation (SD) over the 8 days. The plot above shows the average response for all data points and the ± 2 SD.

Single Cell Performance

Ta-Tagged PBMC




PBMC = peripheral blood mononuclear cell

Ta-Tagged cells were prepared in PBS buffer. The plot above demonstrates replicate measurements from a single sample. Cell transport efficiency will vary depending on cell type and cell stability.


Advantage of the microFAST's Sample Mixing Method

No Mixing Method – Particle Count Stability over 2 hours (n = 20)



Particle count for 50 nm Pt NPs analyzed over a 2 h time period from 20 identically prepared samples in separate vials using the no mixing method. The analysis time was set to ensure the 20 samples took 2 h to complete.


Mixing Method – Particle Count Stability over 2 hours (n = 20)



Particle count for 50 nm Pt NPs analyzed over a 2 h time period from 20 identically prepared samples in separate vials using the mixing method. The analysis time was set to ensure the 20 samples took 2 h to complete.

microFAST SingleCell Features

microFAST SingleCell Features



microFAST SingleCell Autosampler

System

microFAST SingleCell Autosampler

Part Number

MF-SC2-64

Sample Introduction Kit

Single Cell Sample Introduction Kit for Nexion 2000/5000

- CytoNeb 50-64**
CytoNeb 50 Meinhard Glass Nebulizer with PFA Gas Line
- SC-CytoC-64**
CytoSpray Chamber for NexION 2000/5000
- MF-5037-3151-060**
Manual Sampling Line, 0.15 mm ID
- MPP-019-F-PVC**
Orange/Red flared PVC tubing for Micro Peripump (12/pkg)
- ES-2501-PPF2**
CTFE peripump fitting, female, barbed
- S250978**
CytoSpray Chamber Make-Up Gas Line
- MPP-130-PHR**
Gray/Gray Santoprene tubing for Micro Peripump (12/pkg)
- ES-2501-PPM2**
CTFE peripump fitting, male, barbed
- ES-2044-0005**
CytoSpray Chamber Drain Line
- S250871**
PFA CytoNeb Sample Line

ESI Elemental Scientific

Kit

Single Cell Sample Introduction Kit for NexION ICPMS

Part Number

SC-SI-64

*Torch not included with ESI's sample introduction kit. We recommend using PE's fixed 2.0 mm injector UHP quartz torch (N8152428).

FAST Conversion Kit

microFAST SingleCell (MFTM-0370-64)
Total Metals Upgrade Kit

- (MPP-038-F-PVC)**
MP2 Peripump Tubing Flared (org/grn) 12/pkg
- (ICN50-64)**
PFA ICN Nebulizer
- (SC-0318-02)**
Sample Loop, Fluoropolymer, 200µL
- (SC-0318-05)**
Sample Loop, Fluoropolymer, 500µL
- (SC-0318-10)**
Sample Loop, Fluoropolymer, 1 mL
- (SC-5037-4502-C)**
Carrier Probe Line Kit

ESI Elemental Scientific ph: 402.991.7800 | fax: 402.991.7799 | sales@icpms.com | www.icpms.com

Kit

Includes ICN50-64 nebulizer to use with instrument standard spray chamber to run FAST sample analysis on the microFAST SingleCell autosampler

Part Number

MFTM-0370-64



© Elemental Scientific | 7277 World Communications Drive | Omaha, NE 68122
Tel: 1-402-991-7800 | sales@icpms.com | www.icpms.com

B-24222

