

# microFAST SingleCell

**Complete Solution for Single Particle and Single Cell ICPMS Applications**

Automated sample introduction system for Thermo ICPMS




 **Learn More**

**Biological Size Scale**

**Biological Size Scale**

The microFAST SingleCell automation system delivers single nanoparticles or single cells intact to the ICPMS. Once the NPs or cells enter the plasma they are vaporized, atomized, and ionized. The resulting elemental content is detected by the mass spectrometer.



## Single Cell ICPMS


The ability to introduce single cells into an ICPMS and measure the elemental content in each cell, or tagged to each cell, accurately takes a dedicated, well-designed sample introduction system. Having this ability allows for investigators to better understand how much of a specific nanoparticle, metallodrug, or metal-based compound enters the cell. These cells or nanoparticles will vary in size from a few nm's up to a few 100  $\mu$ m's. The typical cell types of interest will vary and with that the stability of the cell-line also varies, such that a gentle, controlled nebulization must be employed in order to not disrupt or lyse the cell.

### Single Cell ICPMS Requirements

- Flexible sample volumes –  $\mu$ L to mL of sample
- Ensure cells stay intact, no cell lysing
- Low pressure sample introduction
- High transport efficiency

Elemental Scientific has developed a complete sample introduction system designed specifically for single cell and nanoparticle applications. This system consists of:

- **microFAST SingleCell Autosampler**
- **CytoNeb** – single cell nebulizer
- **CytoSpray** – single cell spray chamber
- **One-piece Torch** – ICPMS torch for simple, direct connection



## microFAST SingleCell





microFAST SingleCell Automated  
Sample Introduction System for Thermo ICPMS  
Part Number: MF-SC2-73

## CytoNeb and CytoSpray

### CytoNeb


- Meinhard high efficiency nebulizer
- Designed to efficiently nebulize single cells without cell rupturing
- Low internal volume
- Low backpressure (1-50  $\mu\text{L}/\text{min}$  = <50 psi)
- Low dead volume
- Patented inert PFA quick connects for nebulizer gas and samples lines




## CytoNeb and CytoSpray

### CytoSpray

- Spray chamber designed specifically for single cell and nanoparticle applications
- High-transport efficiency
- Separate make-up gas for better transport efficiency
- Includes one-piece ICPMS torch for simple and direct connection to the CytoSpray




## microFAST SingleCell System



Simple Schematic of the setup for measuring Pt in cells


## microFAST SingleCell System



microFAST SingleCell flow path with syringe carrier and sample loading

The microFAST SingleCell system has been built for performance by optimizing the inner diameter (ID) and line lengths to ensure a quick sample transfer from vial to ICP torch.

- Fast sample-to-sample times. For example, at 20  $\mu\text{L}/\text{min}$  flow rate:
  - <3 min, when utilizing a 30 s ICPMS measurement time
  - <4 min, when utilizing a 100 s ICPMS measurement time
- High-flow sample loop washout
- Simple conversion for total metal analysis using FAST system
  - Vacuum or syringe sample loading
  - Micro or large sample volume capabilities



Example of how gold nanoclusters (AuNC) can be used to tag proteins and determine how much of a specific protein is in each cell.

## Thermo Method Setup

### Analysis Options

#### Common

SampleFlow: 10  $\mu\text{L}/\text{min}$

Transport Efficiency:  Auto 28.07 %

Dwell Time: 5.00 ms

#### Advanced

Maximum Probability of Incomplete Measurements: 10.00 %

Maximum Ratio of Filled Slots: 5.00 %

Maximum Duration of Particle Signal: 3.00 ms

Dwell Time Proposal: 30.00 ms

Probability of Incomplete Measurements: 60.00 %

Bin Size relative to Signal Intensity for Dynamic Bin Size: 10.00 %

Auto Background Sigma Multiplier: 5.00

#### Detection Sensitivity

Analyte	Calculate Sensitivity	Sensitivity
195Pt	<input checked="" type="checkbox"/> On	10,979.45 cps/ $\mu\text{g}/\text{L}$

## Thermo Method Setup

### Sample List

Label	Status	Duration (s)	Evaluate	Sample Type	Detection Sensitivity Standard	Transport Efficiency Standard	Rack Number	Vial Numbers
1 BLK1	<input checked="" type="checkbox"/>	30	<input checked="" type="checkbox"/>	UNKNOWN			0	1
2 BLK2	<input checked="" type="checkbox"/>	30	<input checked="" type="checkbox"/>	UNKNOWN			0	1
3 Blank	<input checked="" type="checkbox"/>	30	<input checked="" type="checkbox"/>	DS_BLK			0	1
4 STD1	<input checked="" type="checkbox"/>	30	<input checked="" type="checkbox"/>	DS_STD	New Standard 1		0	2
5 STD2	<input checked="" type="checkbox"/>	30	<input checked="" type="checkbox"/>	TE_STD		New Standard 2	0	3
6 Sample1	<input checked="" type="checkbox"/>	30	<input checked="" type="checkbox"/>	UNKNOWN			1	1
7 Sample2	<input checked="" type="checkbox"/>	30	<input checked="" type="checkbox"/>	UNKNOWN			1	2
8 Sample3	<input checked="" type="checkbox"/>	30	<input checked="" type="checkbox"/>	UNKNOWN			1	3
9 Sample4	<input checked="" type="checkbox"/>	30	<input checked="" type="checkbox"/>	UNKNOWN			1	4
10 Sample5	<input checked="" type="checkbox"/>	30	<input checked="" type="checkbox"/>	UNKNOWN			1	5
11 Sample6	<input checked="" type="checkbox"/>	30	<input checked="" type="checkbox"/>	UNKNOWN			1	6
12 Sample7	<input checked="" type="checkbox"/>	30	<input checked="" type="checkbox"/>	UNKNOWN			1	7
13 Sample8	<input checked="" type="checkbox"/>	30	<input checked="" type="checkbox"/>	UNKNOWN			1	8
14 Sample9	<input checked="" type="checkbox"/>	30	<input checked="" type="checkbox"/>	UNKNOWN			1	9
15 Sample10	<input checked="" type="checkbox"/>	30	<input checked="" type="checkbox"/>	UNKNOWN			1	10

### Define Standards

Standard(s)

Name	Description
10ppb Pt	
50nm Pt	

Transport Efficiency Standard - 50nm Pt

Number per Volume: 400000 particle/mL

Mass: 1.4 fg

Particle Fraction: 195Pt

Description: *Optional*

Standard(s)

Name	Description
10ppb Pt	
50nm Pt	


Elements

Unit: Mass/Volume


Element	Concentration
Pt	10 $\mu\text{g}/\text{L}$

### Example Histogram


Run	Mean Size (nm)
Run 1	45.2
Run 2	44.8
Run 3	45.0
Average	45.0
STD Dev.	0.2
%RSD	0.5



## Single Cell Introduction Kit Performance



Backpressure was recorded for each sample flow rate using DI water, 1X PBS, or 50 mM Tris buffer as the carrier solution. Larger ID tubing can be substituted to achieve lower backpressures.




Backpressure was recorded for varying nebulizer gas flow rates using 10 and 50 µL/min sample flow (DI water as the carrier solution).

## Nanoparticle Performance

### 50 nm Au NPs


Typical Transport Efficiency for 50 nm Au NPs = ~80% or greater



Ten 50 nm Au NPs were prepared under the same conditions and analyzed over an 8-day period. Samples were sonicated before each day's analysis. Data points represent the average response for each sample over the 8 days. Error bars represent  $\pm 1$  standard deviation (SD) over the 8 days. The plot above shows the average response for all data points and the  $\pm 2$  SD.

## Single Cell Performance

### Ta-Tagged PBMC




PBMC = peripheral blood mononuclear cell

Ta-Tagged cells were prepared in PBS buffer. The plot above demonstrates replicate measurements from a single sample. Cell transport efficiency will vary depending on cell type and cell stability.


## Advantage of the microFAST's Sample Mixing Method

### No Mixing Method – Particle Count Stability over 2 hours (n = 20)



Particle count for 50 nm Pt NPs analyzed over a 2 h time period from 20 identically prepared samples in separate vials using the no mixing method. The analysis time was set to ensure the 20 samples took 2 h to complete.

### Mixing Method – Particle Count Stability over 2 hours (n = 20)




Particle count for 50 nm Pt NPs analyzed over a 2 h time period from 20 identically prepared samples in separate vials using the mixing method. The analysis time was set to ensure the 20 samples took 2 h to complete.



**microFAST SingleCell Features**

**microFAST SingleCell Features**



## microFAST SingleCell Autosampler

### System

microFAST SingleCell Autosampler

### Part Number

**MF-SC2-73**

## Sample Introduction Kit

**Single Cell Sample Introduction Kit for Thermo Fisher iCAP Q Series ICP-MS**

**CytoNeb 50-73**  
CytoNeb 50 Meinhard Glass Nebulizer with PFA Gas Line

**T20-73**  
One Piece Torch/Injector for Single Cell 2 mm injector ID, quartz

**SC-CytoC-73**  
CytoSpray Chamber for iCAP Q Series ICP-MS

**ES-2044-0005**  
CytoSpray Chamber Drain Line

**S250827**  
CytoSpray Chamber Make-Up Gas

**ES-2501-PPF2**  
CTFE peripump fitting, female, barbed

**ES-2501-PPM2**  
CTFE peripump fitting, male, barbed

**S250871**  
PFA CytoNeb Sample Line with barbed connection to Peripump

**Q3PT-019-F-PVC**  
PVC Flared Tubing, 3-stop, 0.19 mm ID, (orange/red/orange), 12/pk

**Q3PT-130-PHR**  
Santoprene Tubing, 3-stop, 1.30 mm ID, (gray/gray/gray), 12/pk

**MF-5037-3151-060**  
Manual Sampling Line, 0.15 mm ID

**E**lemental  
**S**cientific

### Kit

Single Cell Sample Introduction Kit for Thermo ICPMS

### Part Number

**SC-SI-73**

## FAST Conversion Kit

**microFAST SingleCell (MFTM-0370-73)**  
**Total Metals Upgrade Kit**

**(MPP-038-F-PVC)**  
MP2 Peripump Tubing Flared (org/grm) 12/pkg

**(ICN50-73)**  
PFA ICN Nebulizer

**(SC-0318-02)**  
Sample Loop, Fluoropolymer, 200µL

**(SC-0318-05)**  
Sample Loop, Fluoropolymer, 500µL

**(SC-0318-10)**  
Sample Loop, Fluoropolymer, 1 mL

**(SC-5037-4502-C)**  
Carrier Probe Line Kit

**E**lemental Scientific ph: 402.991.7800 | fax: 402.991.7799 | sales@icpms.com | www.icpms.com

### Kit

Includes ICN50-73 nebulizer to use with instrument standard spray chamber to run FAST sample analysis on the microFAST SingleCell autosampler

### Part Number

**MFTM-0370-73**



© Elemental Scientific | 7277 World Communications Drive | Omaha, NE 68122  
Tel: 1-402-991-7800 | [sales@icpms.com](mailto:sales@icpms.com) | [www.icpms.com](http://www.icpms.com)

B-24215

