



**Elemental Scientific**

ICP | ICPMS | AA

[www.icpms.com](http://www.icpms.com)

# SAMPLE **TRAX**S5

**Barcoded Labware and Racks**

# SAMPLETRAX S5 PFA Barcoded Labware & Racks

## Barcoded SAMPLETRAX Bottles, Caps and Trays

### Barcoded Caps

#### PFA Vials and Bottles with Barcoded Caps

SAMPLETRAX vials and bottles are available with regular caps or barcoded caps. To order a vial or bottle with the barcoded cap, simply add -BC to the end of the part number.

#### Barcoded Caps - BC

Add -BC to end of vial or bottle part number for barcoded cap  
ex: TXV3-1-BC



- BC

#### Large Volume (>50 mL) PFA Bottles

Volume	OD	Fits Rack
100 mL	48 mm	TXTR-31-100-100-100
250 mL	62 mm	TXTR-31-250-250-250
500 mL	78 mm	TXTR-14-500-500-500

Bottle PN	Barcoded Caps
-----------	---------------

TXV3-100	- BC
TXV3-250	- BC
TXV3-500	- BC



TXV3-250

Add - BC to the end of the Vial Part Number to include barcoded SAMPLETRAX caps

#### Bottle Trays for Large Volume (>50 mL) PFA Bottles

Positions	Bottle OD	Fits Bottle
31	48 mm	TXV3-100
31	62 mm	TXV3-250
14	78 mm	TXV3-500

Tray PN
---------

TXTR-31-100-100-100
TXTR-31-250-250-250
TXTR-14-500-500-500



TXTR-31-250-250-250

\*Note: SAMPLETRAX tray required for large volume PFA bottles

## Barcoded SAMPLETRAX Vials, Bottles and Racks

### Small Volume (<50 mL) PFA Vials and Bottles


Volume	OD	Fits Rack	Vial PN	Barcoded Caps
1 mL	12 mm	TXR-60S-12	TXV3-1	- BC
12 mL	16 mm	TXR-60T-16	TXV3-12	- BC
15 mL	30 mm	TXR-21S-30	TXV3-15	- BC
30 mL	30 mm	TXR-21S-30	TXV3-30	- BC
50 mL	30 mm	TXR-21T-30	TXV3-50	- BC
10 mL	28 mm	TXR-21S-28	TXV3-10	- BC
20 mL	28 mm	TXR-21S-28	TXV3-20	- BC



Add - BC to the end of the Vial Part Number to include barcoded SAMPLETRAX caps

### Racks for Small Volume (<50 mL) PFA Vials and Bottles

Positions	Vial OD	Fits Vial	Rack PN
21	30 mm	TXV3-15, TXV3-30	TXR-21S-30
21	30 mm	TXV3-50	TXR-21T-30
21	28 mm	TXV3-10, TXV3-20	TXR-21S-28
60	12 mm	TXV3-1	TXR-60S-12
60	16 mm	TXV3-12	TXR-60T-16
9	30 mm	TXV3-15, TXV3-30	TXR-ST-2-9-30



Deck PN  
TXD-4

**\*Note: Transparent SAMPLETRAX deck required for small volume PFA vials and bottles**

SAMPLETRAX transparent deck with various racks



**Elemental Scientific**

© Elemental Scientific | 7277 World Communications Drive | Omaha, NE 68122  
Tel: 402-991-7800 | [sales@icpms.com](mailto:sales@icpms.com) | [www.icpms.com](http://www.icpms.com)