


External Location

- 1.7 m (5.5') PTFE-shielded control cable
- Allows regulator to be located outside hood to eliminate corrosion
- 110V or 208-240V supply

Control

- +2°C control up to 250°C maximum
- Multiple alarms for maximum safety
- Non-programmable with timer shutoff
- Programmable version with multiple temperature steps



Single Controller Regulator

- Provides heating control for one ANALAB device
- Temperature readout from device

Dual Controller Regulator

- Provides heating control for one ANALAB device
- Temperature readout from device
- Temperature readout from external probe or control from external probe