

AutoAlign Arm

- Self-realigns after encountering capped sample tube
- Metal-free fluoropolymer probe aligner

Sample Containers

- <1 mL up to 500 mL

PTFE Autosampler Deck

- Chemically resistant for long lifetime

Position Autocorrection

- Adapts to every position automatically
- Automatically retries position if there is an obstacle
- Keeps analysis running, after years of use

Parabolic Probe

- Robust, heavy-walled, self-supporting
- Upgradable to PFA Fluorall Probe

Fluoronic Drive

- True fluoropolymer mechanical action
 - Magnetically coupled linear drive
- Chemically resistant and long-lived

Dual Overflowing Rinse Stations

- Ultra-low cross-contamination
- Independently pumped for low rinse consumption

SampleSense *FAST* & prep*FAST* Option

- Automated, ultra-high sample throughput with SampleSense *FAST*
- Autocalibration, autodilution, and automated sample loading with SampleSense prep*FAST*

