

Septum-Piercing Autosampler

Highest Speed of Operation

- Higher speed analysis
- Keeps analysis running, even after obstruction

Septum-Piercing Probe

- High purity sapphire probe
- For piercing PTFE, silicon, Viton® & butyl septa

PTFE Autosampler Deck

- Chemically resistant for long lifetime

Position Autocorrection

- Adapts to every position automatically
- Automatically retires position if there is an obstacle

Sample Containers

- Septa sealed and unsealed vessels
- < 1mL up to 500 mL

Dual Overflowing Rinse Stations

- Ultra-low cross-contamination
- Optional organic solvent rinse
- Independently pumped for low rinse consumption

SampleSense FAST & prepFAST Ready

- Automated ultra-high sample throughput with SampleSense FAST
- Autocalibration, autodilution and automated sample loading with SampleSense prepFAST

