

# TRAXOil Features Diagram

## TRAXGrip

- Purpose-built for consistent, repeatable and reliable oil sample handling through workflow

## TRAXDisplay

- Powered by ESI's LabSymphony Software with direct LIMS, analytical software export and unlimited method storage
- Integrated touchscreen control for method creation, workflow execution, and system monitoring

## TRAXSpin

- Automated vortex mixing keeps additives, particulates, and wear metals evenly suspended for representative sampling.
- Improves aliquot consistency and analytical reproducibility from sample-to-sample

## LinkGauge

- Infrared sensing confirms cap removal and process completion in real time
- Dynamic liquid-level detection improves aspiration precision and minimizes sample carryover

## TRAXLink

- Uncaps and recaps samples
- Verifies sample identity with integrated barcode scanning

## Heating Block

- Uniformly warms samples to reduce viscosity, improving aliquot accuracy and aspiration across different oil sample types

## Waste Management Cabinet

- Simplifies waste line management
- Leak sensors alert users of spills and overfills

