

TRAXStation Features

ULPAClean 10 Filters

- Removes 99.9995% of airborne particles

TRAXDualDisplay

- Easily run and build methods with ESI's innovative software

TRAXShield

- Clean sample environment
- Protects lab from harmful fumes
- Easily access samples from all four sides

TRAXGripper

- Transfers sample tubes between source rack and TRAXLink

Waste Management Cabinet

- Clean exchange of waste bottle – no drips
- Simplified waste line management
- Level and leak sensors

StationPump

- Post-mounted high pressure syringe module for easy syringe visibility
- Precise aliquots as low as 10 μ L
- Rapid and precise dispensing

LuerProbe

- Attaches to and dispenses through filters

Dual Rinse Station

- Ultra-low cross-contamination
- Low rinse consumption

Filter Ejection Station

- Automatically ejects used filters

TRAXLink

- Automatically uncaps sample vials, scans barcodes, and re-caps

Filter Disposal Bin

- Collects used filters for later disposal

