

Optional ULPAclean 10 Filter

- Removes 99.9995% of airborne particles

AutoAlign Arm

- Self-realigns after encountering capped sample tube

Exhausted Enclosures

- Three side easy access doors

Carrier and Working Solution

Sample Containers

- 0.5 mL vials to 30 mL tubes
- Microplates

Syringes

- High precision syringes
- Stable injection at low flow rates (e.g. 10 μ L/min)

PTFE Autosampler Deck

CytoSpray

- Specifically designed for single cell and nanoparticle applications
- Designed for high-transport efficiency

One-piece ICPMS Torch

- Simple and direct connection to CytoSpray

CytoNeb

- Low backpressure and dead volume
- Unmatched nebulization efficiency

Magnetic SnapValves

Fluoronic Rail

- True fluoropolymer mechanical action
- Magnetically coupled linear drive
- Chemically resistant and long-lived

Tapered Probe

Optional FluoroWand Arm

Dual Pumped Rinse Station

- Low rinse consumption
- Two rinse solutions

