



### Dual Calibrated Sample Loops

- Software controlled sample volumes of 5  $\mu$ L to 1000  $\mu$ L
- Inject sample while next sample loads into other loop
- Magnetic coupling technology

### Magnetic Snap Valves

### Exhausted Enclosures

- Three side easy access doors

### Carrier and Working Solution

### FluoroWand Racks

- Custom racks for FluoroWand compatible vials, bottles and jars
- Elevated to locate vessel tops at consistent height for FluoroWand operation

### Syringes

- High precision syringes
- Stable injection at low flow rates (e.g. 10  $\mu$ L/min)

### Optional ULPAclean 10 Filter

Removes 99.9995% of airborne particles

### Fluoronetic Rail

- True fluoropolymer mechanical action
- Magnetically coupled linear drive
- Chemically resistant and long-lived

### AutoAlign Arm

- Self-realigns after encountering capped sample tube

### FluoroWand

- Removes and replaces sample covers with a vacuum suction cup
- Minimizes risk of contamination
- Prevents sample evaporation and gas phase reactions of samples

### Dual Pumped Rinse Station

- Low rinse consumption
- Two rinse solutions

### PFA Sample Containers

- <1 mL up to 100 mL

### PTFE Autosampler Deck