


CytoSpray

- Specifically designed for single cell and nanoparticle applications
- Designed for high-transport efficiency



CytoNeb

- Low backpressure and dead volume
- Unmatched nebulization efficiency

Carrier and Working Solution

Sample Containers

- 0.5 mL vials to 30 mL tubes
- Microplates

Syringes

- High precision syringes
- Stable at low flow rates (e.g. 10 $\mu\text{L}/\text{min}$)

Optional ULPA Filter or Air Shower

Exhausted Enclosure

- Three side easy access doors

Fluoronetic Z-rail

- True fluoropolymer mechanical action
 - Magnetically coupled linear drive
 - Chemically resistant and long-lived

AutoAlign Arm

- Self-realigns after encountering capped sample tube

Tapered Probe

Dual Rinse Station

- Low rinse consumption

PTFE Autosampler Deck

