

prepFAST Station Features

Dual Lumen Probe

- Two channels integrated in one probe
- One channel for picking up sample aliquot
- Second channel for mixing diluted sample in the destination vial
- Second probe can also be used to stir up sample prior to aliquot

Sample Containers

- Hundreds of compatible racks and vials
- Compatible with microplates

Sample Syringe

- Used to pick up sample aliquot from source vial

Integrated Touchscreen

- Easy-to-use software
- Stores an unlimited number of custom methods

AutoAlign Arm

- Self-aligns after encountering capped sample tube
- Metal-free fluoropolymer probe aligner

Dual Overflowing Rinse Stations

- Ultra-low cross-contamination
- Independently pumped for low rinse consumption

Diluent Syringe

- Used to add diluent in Destination vial

Mass Flow Controller

- Used to mix samples after each sample is prepared
- Can be used to mix samples prior to aliquot to ensure homogenous sample uptake

