

prepFAST X

AutoAlign Arm

- Self-realigns after encountering capped sample tube
- Metal-free fluoropolymer probe aligner

Sample Containers

- Flexible rack configurations
- Samples from micro vials to large bottles

Autocorrecting x-y-z

- Automatically goes to the correct position every time
- Automatically retries position if an obstacle is encountered

Highest Speed of Operation

- Corrects for stall events, even after years of use
- Keep analysis running, even after obstruction

Inline Syringe Sample Preparation

- Sample dilution
- Internal standard addition
- Overage dilution

Large Volume Reagent Bottles

Parabolic Probe

- Robust, heavy-walled, self-supporting

Optional Vented Enclosure

- ULPA filter option

Dual Overflowing Rinse Stations

- Ultra-low cross-contamination
- Optional organic solvent rinse
 - Independently pumped for low rinse consumption

Vacuum Pump

- Vacuum loaded sampling
- Rapid rinsing & cleaning

Chemically Resistant

- Designed for exceptionally long lifetime
- No metallic components above samples

SampleSense Valve Module

- Automatic optically-sensed sample loading
- Single valve 400x dilution range

