

seaFAST S2

Electronics in Rear Plenum

Exhausted Enclosure

Built-in to provide maximum contamination control

Sample Containers

<1 mL up to 500 mL

Over-Flowing Dual Rinse Station

Minimize sample contamination

Ultraclean Valve System

Completely enclosed liquid flow path

Syringe-Driven Chemistry

Inline or offline, ultrapure preconcentration and matrix removal. Direct mode with static inline dilution.

Bookshelf Reagent Bottles

Maximizes storage for all reagents

