

## 8DXCi Metal-free AutoBench

**All-in-one Solution For Your Lab**  
Compact Analytical Bench With  
Integrated 8-rack Autosampler



### Features & Benefits

- **All-in-one solution**
  - Extend analysis runs and maximize efficiency with 8-rack autosampler
  - Improve lab safety with innovative waste containment cabinet
  - Seamlessly integrates with ICP or ICPMS
- **Metal-free components**
  - Minimize risk of contamination typically caused by metal components
  - Enhance sensitivity, improve detection limits, and improve reproducibility
  - Extend system lifespan and reduce maintenance costs
- **Compact footprint**
  - Save valuable lab floorspace with compact system footprint that fits an 8-rack autosampler, computer, ICP/ICPMS, waste bottle, vacuum pump, and more
  - System Dimensions:

**Width:** 160 cm (63 in)

**Height:**

With Monitor - 195 cm (76.5 in)

Floor to Work Surface - 80 cm (31.5 in)

*Height will vary based on instrument*


**Depth:** 86 cm (34 in)

*Also available in 4-rack configuration*

# Features Diagram

## Dual Monitor Display

Simultaneously monitor your ICPMS and analyze incoming datapoints




## Clean3V Enclosure

Vertical sliding doors provide convenient access to samples from three sides


## ULPAClean 10

Removes 99.9995% of particles from enclosure



## Consumables Drawer

Convenient additional storage for cones, nebulizers, & more



## Vacuum Pump Silencing Cabinet

Reduces vacuum pump noise and vibration

## Waste Management Cabinet

Eliminates hassle with waste lines and prevents large messes from spills or overfills

## Storage Cabinet

Provides convenient storage for syringes, valves, bottles, and prepFAST automation accessories

