

# IsotopeBench



## COMPACT METAL-FREE LAB BENCH FOR ISOTOPE2 & APEX OMEGA

The IsotopeBench is an all-in-one solution for your laboratory. IsotopeBench is built from metal-free materials and comfortably holds the Isotope2 and Apex Omega, with additional storage cabinets and drawers.

Engineered with a compact footprint to optimize space utilization, this system enhances laboratory efficiency while offering ample storage capacity. It's innovative waste and exhaust management system is vital to improving the safety of labs.

### Key Benefits:

- Designed for Isotope2 & Apex Omega systems
- Metal-free and corrosion-free
- Improve lab safety with waste & exhaust management
- Welded polypropylene

<b>Part Number:</b> AB-Isotope2
<b>Width:</b> 162 cm (64 in)
<b>Height:</b> 122 cm (48 in)
<b>Depth:</b> 61 cm (24 in)
<b>Weight:</b> 246 kg (581 lbs)

*Easily mobile for exchanging with other sample intro systems such as ESL laser ablation systems.*



Learn More



© Elemental Scientific | 7277 World Communications Drive | Omaha, NE 68122  
Tel: 1-402-991-7800 | sales@icpms.com | www.icpms.com

F-24109

