


Cold P'Tips Cleaner

*Cold conditioning device, simple and effective
for your pipette tips*



Filling position

- Don't take the risk of contaminating your samples and/or reagent you collect by pipetting.
- Pre-conditioning the inside and outside of your pipette tips in a single effective operation, reproducible, economic.
- Possibility to clean any type of pipette tips.



Cleaning Holder

96 pipettes tips transferred at once in the cleaning holder.

Device made in fluoropolymer materials (PTFE, PFA) can be used with a large number of reagents.



Drain position